

An ISO 13485:2016 and
ISO 9001:2015 Certified Company



GI Endoscopy

Company Profile

Endo-Med Technologies Pvt. Ltd. is having its corporate office at New Delhi, capital city of India and manufacturing facility at Noida (National Capital Region) engaged in development, manufacturing and globally distribution of minimal invasive medical devices in the field of Gastro Endoscopy, Laparoscopy, and Bronchoscopy. Backed by more than 25 years of experience. Endo-Med is committed to serve the customers with quality products at affordable price and service to the utmost satisfaction of the customers.

Quality Policy

The members of "team Endo-Med" are experienced with the largest companies across the world and exposed to trans-continent markets, thus possessing the highest level of expertise in their respective fields. The company is an ISO 13485:2016 certified organization which is deeply involved in the manufacturing of the best quality medical products under the best quality management systems. With guidance from the renowned Doctors in India and abroad associated with Endo-Med, we develop the product keeping in mind, the every minute expectation of the user from the product thus making our product among the most user friendly for a medical practitioner in their procedural room. All Endo-Med Products are "CDSCO" Certified and are manufactured as per the guidelines of GMP*.

Philosophy

At Endo-Med, Quality is the way of life and operations. Stage wise rigorous inspection is carried out right from procurement of raw material to delivery of the finished product. All of our employees are committed to constantly improve the quality management system.


The company believes in long term relationship and mutual growth as objective. At Endo-Med, quality is a culture and service a tradition.



FOREIGN BODY NET BASKETS SINGLE USE

Product Features

- Rotatable
- Wide Opening
- Extra Strong Net
- Smooth Operation Even in Loops
- Easy Grip Handle




Cat No.	Description	Sheath OD	Length
EFB-0730230	Net Type 30 mm	2.3 mm	230 cm
EFB-0740230	Net Type 40 mm	2.3 mm	230 cm


POLYPECTOMY SNARES

Product Features

- Braided Wire
- Adjustable Size
- Easily Rotatable
- 360° Rotatable



Cat No.	Type	Snare Size	Sheath Size
EPS- 0125230	Oval Shape	25 mm	2.4 mm x 230 cm
EPS- 0130230	Oval Shape	30 mm	2.4 mm x 230 cm
EBS - 0140230	Oval Shape	40 mm	2.4 mm x 230 cm
EBS - 0150230	Oval Shape	50 mm	2.4 mm x 230 cm
EPS- 0230230	Hexagonal Shape	30 mm	2.4 mm x 230 cm



Endo-Med Technologies Pvt. Ltd.



Corporate Off. : Off. No. 607, 6th Floor, DDA Building,
Plot No. 4, Laxmi Nagar Distt. Centre, Delhi-110092. INDIA

Factory : A-138, Sector 63, Noida - 201301 U.P.

Phone: +91-11-22048222, 42448449

Customer Care No. :  +91-7840070477

E-mail: sales@endomed.in, **Web** : www.endomed.in

 Endo-Med Technologies Pvt. Ltd.  endomed_technologies

Distributed By :