


Gastrostomy Feeding Tube

Product Features

- Medical Grade 100% Silicon Tube
- Thin Walls Tube
- Large Inner Lumen
- Non-traumatic Distal Tip



| Cat No. | Tube OD | Tube Length | Balloon Volume |
|-------------|---------|-------------|----------------|
| EPG-2112024 | 12 Fr | 24 cm | 5 ml |
| EPG-2114024 | 14 Fr | 24 cm | 5 ml |
| EPG-2116024 | 16 Fr | 24 cm | 20 ml |
| EPG-2118024 | 18 Fr | 24 cm | 20 ml |
| EPG-2120024 | 20 Fr | 24 cm | 20 ml |
| EPG-2122024 | 22 Fr | 24 cm | 20 ml |
| EPG-2124024 | 24 Fr | 24 cm | 20 ml |