


MSD - Multi Stage Dilation

DILATION BALLOON

Esophageal, Pyloric / Colonic, Biliary, Achalasia

Product Features

- 3 stage Dilatation in one Balloon
- Transparent Material
- Metal Marks For Easy Positioning
- Guide Wire Compatible
- Non-Compliant Balloon
- Fast Inflation & Deflation



Cat No.	Application	Balloon OD	Balloon Length	Total Length	Channel Id
EBD-M106230	Biliary/Esophagus/ Pylorus/Colon	6-7-8 mm	55 mm	230 cm	>2.8 mm
EBD-M108230	Biliary/Esophagus/ Pylorus/Colon	8-9-10 mm	55 mm	230 cm	>2.8 mm
EBD-M110230	Esophagus/Pylorus/ Colon	10-11-12 mm	80 mm	230 cm	>3.2 mm
EBD-M112230	Esophagus/Pylorus/ Colon	12-13.5-15 mm	80 mm	230 cm	>3.2 mm
EBD-M115230	Esophagus/Pylorus/ Colon	15-16.5-18 mm	80 mm	230 cm	>3.2 mm
EBD-M118230	Esophagus/Pylorus/ Colon	18-19-20 mm	80 mm	230 cm	>3.2 mm
EBD-0430080*	Achalasia	30 mm	80 mm	90 cm	N/A
EBD-0435080*	Achalasia	35 mm	80 mm	90 cm	N/A
EBD-0440080*	Achalasia	40 mm	80 mm	90 cm	N/A

*Non- Multi Stage Balloons