


V-JECTOR[®]

ENDOSCOPY SCLERO THERAPY NEEDLE

Product Features

- Lock Mechanism
- Easy Grip Handle
- 45° Beveled Stainless Steel Needle
- Metal Tip
- Smooth operation even at difficult locations
- PTFE inner & outer Tubes



Cat No.	Needle Gauge	Needle Length	Sheath OD	Sheath Length
ESN-0121180	21G	5 mm	2.3 mm	180 cm
ESN-0123180	23G	5 mm	2.3 mm	180 cm
ESN-0121230	21G	5 mm	2.3 mm	230 cm
ESN-0123230	23G	5 mm	2.3 mm	230 cm