

**Vu-Max**<sup>®</sup>  
6+1 BANDS LIGATOR

**MULTI BAND LIGATOR**

**Product Features**

- Maximum Endoscopic View
- Red Color Trail Band
- Penultimate Blue Color Band
- Six Usable Bands
- Easy Irrigation / Flushing Port
- Click Sound With Every Band Release



Cat No.	Scope OD	Band per Barrel	Inclusion
EBL-0106085-7 (Full set)	9 mm to 11.5 mm	Six + One (Red)	Handle+ Barrel+Loading Catheter + Flushing Tube
EBL-0206085-7 ( Mini Pack)	9 mm to 11.5 mm	Six + One (Red)	Barrel+Loading Catheter