


EUS & EBUS FNA SYSTEM

Product Features

- Superior Needle Visibility
- Sharp Needle Atraumatic Puncture
- Large Needle opening profile ensure better specimen



Cat No.	Needle G	Needle Length	OD	Length	For Scope
EUS-0119142	19 G	0-80 mm adjustable	1.8 mm	142 cm	Olympus/Pentax/Fuji EUS
EUS-0122142	22 G	0-80 mm adjustable	1.6 mm	142 cm	Olympus/Pentax/Fuji EUS
EUS-0222073	22 G	0-50 mm adjustable	1.6 mm	73 cm	Olympus EBUS
EUS-0322073	22 G	0-50 mm adjustable	1.6 mm	73 cm	Pentax EBUS