


EXPANDER[®]

ESOPHAGEAL DILATOR SET

Product Features

- Set of Eight & Sixteen Dilators
- Total Length 85 cm
- Marking Every 1 cm
- Built-in Radiopaque Metal Ring
- Highly Smooth and Flexible



Cat No.	Set Includes	Dilator Sizes
EDS-0108085	Eight Dilators + One Spring Tip Wire	5,7,9,11,13,15, 17,19 mm x 85 cm
EDS-0116085	Sixteen Dilators + Spring Tip Wire	5,6,7,8,9,10,11,12,13,14,15,16,17,18, 19,20 mm x 85 cm