


Product Features

- Nitinol Stent Highly Radiopaque
- Covered From Inside
- Easy Retrieval Using Drawstring
- Dynamic Radial Force for Maximum Lumen Size
- Dumbbell Shape for Anti-Migration
- Flexible Construction to Conform to Anatomy
- X-Ray Visible Marks
- Scaled Marking on Delivery System



Cat No.	Stent OD	Flare OD	Stent Length	Delivery Size
EMS-0118008	18 mm	22 mm	8 cm	7 mm x 75 cm
EMS-0118010	18 mm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012	18 mm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015	18 mm	22 mm	15 cm	7 mm x 75 cm
EMS-0118010(AR)	18 mm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012(AR)	18 mm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015(AR)	18 mm	22 mm	15 cm	7 mm x 75 cm

*Other sizes are available on request