



PALLIATOR[®] Pseudocyst Stent

Product Features

- Nitinol Stent
- Fully Covered
- With Lasso
- Braided Tubing of Delivery System
- Three Marking on Delivery System for Reliable Positioning



Cat No.	Stent Type	Stent OD	Stent Length	Delivery Size
EMS-0414020	Covered	14 mm	2 cm	10Fr x 180cm
EMS-0414030	Covered	14 mm	3 cm	10Fr x 180cm
EMS-0414040	Covered	14 mm	4 cm	10Fr x 180cm
EMS-0416020	Covered	16 mm	2 cm	10Fr x 180cm
EMS-0416030	Covered	16 mm	3 cm	10Fr x 180cm
EMS-0416040	Covered	16 mm	4 cm	10Fr x 180cm