



Endo-Med[®]

A Resurgence to life

Gastroenterology


Endoscopy

Rota-Net[™]

FOREIGN BODY NET BASKETS SINGLE USE

Product Features

- Rotatable
- Wide Opening
- Extra Strong Net
- Smooth Operation Even in Loops
- Easy Grip Handle



Cat No.	Description	Sheath OD	Length
EFB-0730230	Net Type 30 mm	2.3 mm	230 cm
EFB-0740230	Net Type 40 mm	2.3 mm	230 cm