


## RETRIEVAL BAG

### Product Features

- The bag body of high strength with good biocompatibility is TPU material, Transparent soft bag body
- There are a variety of specifications bag body to choose from
- Meet the demand of different operation
- When the fetch, the automatic tightening of mouth, to prevent overflow
- The operation is simple, The efficiency of safety
- Can prevent cancer to grow and tissue fluid permeability caused by inflammation, Reduce the risk of tissue damage




### Specifications

RE10-250 $\Phi$ 10×250ml	RE10-M $\Phi$ 10×200ml
RE10-350 F 10×350ml	RE10-L F 10×300ml
RE10-800 F 10×800ml	

## SUCTION & IRRIGATION SET

### Product Features

- Ergonomically designed handpiece can be easily operated by either hand
- Unblocked channel, good flowing performance
- Stainless steel probe with position function, easy to fix and draw out
- Flip designed rear port for easy access inside chamber, can be operated by one hand, quick close and reduce air leakage
- Transparent headpiece, unobstructed visualization
- Resistant structure and reliable quality
- Clear channel, consistent diameter with self-flush function
- Smoke evacuation and filtration, smoke flow rate is adjustable



### Specifications

SI06-360T $\Phi$ 5×360mm
SI06-430T F 5×430mm