


EBUS FNA SYSTEM

Product Features

- Superior Needle Visibility
- Sharp Needle Atraumatic Puncture
- Large Needle opening profile ensure better specimen



Cat No.	Needle G	Needle Length	OD	Length	For Scope
EUS-0222073	22 G	0-50 mm adjustable	1.6 mm	73 cm	Olympus EBUS
EUS-0322073	22 G	0-50 mm adjustable	1.6 mm	73 cm	Pentax EBUS


*Other sizes are available on request

DILATION BALLOON

AIRWAY DILATION BALLOON

Product Features

- Transparent Material
- Metal Marks For Easy Positioning
- Guide Wire Compatible
- Non-Compliant Balloon
- Fast Inflation & Deflation



Cat No.	Application	Balloon OD	Balloon Length	Total Length	Channel Id
EBD-0106120	Airway	6 mm	4 cm	120 cm	> 2.8 mm
EBD-0108120	Airway	8 mm	4 cm	120 cm	> 2.8 mm
EBD-0110120	Airway	10 mm	4 cm	120 cm	> 2.8 mm
EBD-0112120	Airway	12 mm	4 cm	120 cm	> 2.8 mm

*Other sizes are available on request