


FOREIGN BODY FORCEPS RE-USABLE

Product Features

- Wide Opening
- Extra Strong Grip
- Grinded Distal Spiral End
- Auto-clavable
- Smooth Operation Even in Loops
- Easy Grip Handle



Cat No.	Description	OD	Length
EFB-R123180	Rat Tooth Type	2.3 mm	180 cm
EFB-R123230	Rat Tooth Type	2.3 mm	230 cm
EFB-R223180	Alligator Type	2.3 mm	180 cm
EFB-R223230	Alligator Type	2.3 mm	230 cm
EFB-R323180	Rat Tooth Alligator Type	2.3 mm	180 cm
EFB-R323230	Rat Tooth Alligator Type	2.3 mm	230 cm

*Other sizes and configurations are available on request