

GI Endoscopy

360° Rotatable



GeM
Government
e Marketplace



ISO 13485-2016



Company Profile

Endo-Med Technologies Pvt. Ltd. is headquartered in New Delhi, India, with its manufacturing facility located in Noida, part of the National Capital Region. We specialize in the development, manufacturing, and global distribution of minimally invasive medical devices, focusing on Gastro Endoscopy, Laparoscopy, and Bronchoscopy. With over 27 years of experience, Endo-Med is dedicated to providing high-quality products at affordable prices while ensuring exceptional customer satisfaction through reliable service.

Quality Policy


At Endo-Med, our team consists of highly skilled professionals with extensive experience in global markets, bringing expertise from some of the largest companies in the industry. As an ISO 13485:2016 certified organization, we adhere to rigorous quality management systems in the manufacturing of top-tier medical products. Guided by renowned doctors both in India and internationally, we design products that meet the specific, detailed needs of medical practitioners, ensuring ease of use in clinical environments. All Endo-Med products are CDSCO certified and manufactured in accordance with GMP guidelines.

Philosophy

Quality is at the heart of everything we do at Endo-Med. We conduct rigorous, stage-wise inspections from raw material procurement to the final delivery of products, ensuring that every aspect meets the highest standards. Our employees are committed to continuously improving our quality management system.

We believe in fostering long-term relationships with our customers, focusing on mutual growth and success.

Endo-Med Manufacturing Facility at Noida



At Endo-Med,

Quality is more than a process—It's a culture;

Service is not just a function, but a tradition.

Smart-BiteTM

BIOPSY FORCEPS SINGLE USE

Product Features

- Wide Opening
- Large Cup Size
- Cups with Back Eye
- With PTFE Sheath covering
- Sterile - Single Use
- Smooth Operation Even in Loops
- Easy Grip Handle




Cat No.	Description	OD	Length
EBF-S110290	Serrated Cup	1.0 mm	290 cm
EBF-S118180	Serrated Cup	1.8 mm	180 cm
EBF-S123180	Serrated Cup	2.3 mm	180 cm
EBF-S123230	Serrated Cup	2.3 mm	230 cm
EBF-S223180	Serrated Cup with needle	2.3 mm	180 cm
EBF-S223230	Serrated Cup with needle	2.3 mm	230 cm
EBF-S123180-H	Oval Cup Hot Biopsy	2.3 mm	180 cm
EBF-S123230-H	Oval Cup Hot Biopsy	2.3 mm	230 cm

*Other sizes and configurations are available on request

FOREIGN BODY FORCEPS RE-USABLE

Product Features

- Wide Opening
- Extra Strong Grip
- Grinded Distal Spiral End
- Auto-clavable
- Smooth Operation Even in Loops
- Easy Grip Handle



Cat No.	Description	OD	Length
EFB-R123180	Rat Tooth Type	2.3 mm	180 cm
EFB-R123230	Rat Tooth Type	2.3 mm	230 cm
EFB-R223180	Alligator Type	2.3 mm	180 cm
EFB-R223230	Alligator Type	2.3 mm	230 cm
EFB-R323180	Rat Tooth Alligator Type	2.3 mm	180 cm
EFB-R323230	Rat Tooth Alligator Type	2.3 mm	230 cm

*Other sizes and configurations are available on request



Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy

GI Endoscopy
GI Endoscopy


GI Endoscopy



FOREIGN BODY NET BASKETS SINGLE USE

Product Features

- 360° Rotatable
- Wide Opening
- Extra Strong Net
- Smooth Operation Even in Loops
- Easy Grip Handle




Cat No.	Description	Sheath OD	Length
EFB-0730230	Net Type 30 mm	2.3 mm	230 cm
EFB-0740230	Net Type 40 mm	2.3 mm	230 cm

POLYPECTOMY SNARES

Product Features

- Braided Wire
- Adjustable Size
- Easily Rotatable
- 360° Rotatable





Cat No.	Type	Snare Size	Sheath Size
EPS- 0125230	Oval Shape	25 mm	2.4 mm x 230 cm
EPS- 0130230	Oval Shape	30 mm	2.4 mm x 230 cm
EBS - 0140230	Oval Shape	40 mm	2.4 mm x 230 cm
EBS - 0150230	Oval Shape	50 mm	2.4 mm x 230 cm

G-clip[™]

HEMOCLIP SINGLE USE

Product Features

- Pre-loaded clip
- 360° Rotatable
- Wide Opening
- Repeated open & close possible
- Smooth Operation at difficult location
- Easy Grip Handle




Cat No.	Description	Opening	OD	Length
EHS-0123230	Pre-loaded	11 mm	2.3 mm	230 cm
EHS-0223230	Pre-loaded	13 mm	2.3 mm	230 cm
EHS-0323230	Pre-loaded	16 mm	2.3 mm	230 cm



ENDOSCOPY SCLERO THERAPY NEEDLE

Product Features

- Lock Mechanism
- Easy Grip Handle
- 45° Beveled Stainless Steel Needle
- Metal Tip
- Smooth operation even at difficult locations
- PTFE inner & outer Tubes



Cat No.	Needle Gauge	Needle Length	Sheath OD	Sheath Length
ESN-0121180	21G	5 mm	2.3 mm	180 cm
ESN-0123180	23G	5 mm	2.3 mm	180 cm
ESN-0121230	21G	5 mm	2.3 mm	230 cm
ESN-0123230	23G	5 mm	2.3 mm	230 cm



Multi Band Ligator (7 Bands)

Product Features

- Maximum Endoscopic View
- Red Color Trail Band
- Penultimate Blue Color Band
- Six Usable Bands
- Easy Irrigation / Flushing Port
- Click Sound With Every Band Release




Cat No.	Scope OD	Band per Barrel	Inclusion
EBL-0106085-7 (Full set)	9 mm to 11.5 mm	Six + One (Red)	Handle+ Barrel+Loading Catheter + Flushing Tube
EBL-0206085-7 (Mini Pack)	9 mm to 11.5 mm	Six + One (Red)	Barrel+Loading Catheter

Vari-fix[®]

MULTI BAND LIGATOR (6 BANDS)

Product Features

- Six Usable Bands
- Different Color Penultimate Band
- Easy Irrigation / Flushing Port



Cat No.	Scope OD	Band per Barrel	Inclusion
EBL-0106085-6 (Full set)	9 mm to 11.5 mm	Six	Handle+ Barrel+Loading Catheter + Flushing Tube
EBL-0206085-6 (Mini set)	9 mm to 11.5 mm	Six	Barrel+Loading Catheter

Vari-fix[®] Ultra

MULTI BAND LIGATOR (7 BANDS)

Product Features

- Seven Usable Bands
- Extremely Fine Shaped Bands
- Different Color Penultimate Band
- Easy Irrigation / Flushing Port
- Micro Beads Based Thread
- Very Thin Natural Yarn



Cat No.	Scope OD	Band per Barrel	Inclusion
EBL-0107085 (Full Set)	9 mm to 11.5 mm	Seven	Handle+ Barrel+Loading Catheter + Flushing Tube
EBL-0207085 (Mini Set)	9 mm to 11.5 mm	Seven	Barrel+Loading Catheter




DILATION BALLOON

Esophageal, Pyloric, Colonic, Biliary & Achalasia

Product Features

- 3 stage Dilatation in one Balloon
- Transparent Material
- Metal Marks For Easy Positioning
- Guide Wire Compatible
- Non-Compliant Balloon
- Fast Inflation & Deflation



Catalogue Nos.	Application	Balloon OD	Balloon Length	Total Length	Minimum Scope Channel	With / Without GW
EBD-M106230	Biliary/Esophagus/Pylorus/Colon	6-7-8mm	80mm	230cm	>2.8mm	Without GW
EBD-M108230	Biliary/Esophagus/Pylorus/Colon	8-9-10mm	80mm	230cm	>2.8mm	Without GW
EBD-M110230	Biliary/Esophagus/Pylorus/Colon	10-11-12mm	80mm	230cm	>2.8mm	Without GW
EBD-M112230	Biliary/Esophagus/Pylorus/Colon	12-13.5-15mm	80mm	230cm	>2.8mm	Without GW
EBD-M115230	Biliary/Esophagus/Pylorus/Colon	15-16.5-18mm	80mm	230cm	>3.2mm	Without GW
EBD-M118230	Biliary/Esophagus/Pylorus/Colon	18-19-20mm	80mm	230cm	>3.2mm	Without GW
EBD-M206230	Biliary/Esophagus/Pylorus/Colon	6-7-8mm	55mm	230cm	>2.8mm	With GW
EBD-M208230	Biliary/Esophagus/Pylorus/Colon	8-9-10mm	55mm	230cm	>2.8mm	With GW
EBD-M210230	Biliary/Esophagus/Pylorus/Colon	10-11-12mm	55mm	230cm	>2.8mm	With GW
EBD-M212230	Biliary/Esophagus/Pylorus/Colon	12-13.5-15mm	55mm	230cm	>2.8mm	With GW
EBD-M215230	Biliary/Esophagus/Pylorus/Colon	15-16.5-18mm	55mm	230cm	>3.2mm	With GW
EBD-M218230	Biliary/Esophagus/Pylorus/Colon	18-19-20mm	55mm	230cm	>3.2mm	With GW
*EBD-0430080	Achalasia	30mm	80mm	90cm	N/A	Without GW
*EBD-0435080	Achalasia	35mm	80mm	90cm	N/A	Without GW
*EBD-0440080	Achalasia	40mm	80mm	90cm	N/A	Without GW

*Single Stage Balloons




EUS & FNA / FNB SYSTEM




Product Features

- Superior Needle Visibility
- Sharp Needle Atraumatic Puncture
- Large Needle opening profile ensure better specimen

FNA



FNB




Catalogue Nos.	Needle G	Max. Needle Length	Sheeth OD	Sheeth Length	FNA / FNB	Scope Compatability
EUS-0119140	19G	80mm	1.8mm	140cm	FNA	Olympus, Pentax, Fuji
EUS-0122140	22G	80mm	1.6mm	140cm	FNA	Olympus, Pentax, Fuji
EUS-0219140	19G	80mm	1.8mm	140cm	FNB	Olympus, Pentax, Fuji
EUS-0222140	22G	80mm	1.6mm	140cm	FNB	Olympus, Pentax, Fuji
EUS-0225140	25G	80mm	1.6mm	140cm	FNB	Olympus, Pentax, Fuji

AUTOMATIC BIOPSY NEEDLE

BIOPSY GUN

Product Features

- Trucut Type Biopsy Sample
- Automatic Biopsy Sampling
- Marking Every 1 CM
- Easy Gripe Handle




Cat No.	Biopsy Size	OD	Length
EBN-0114010	20 mm x 0.8 mm	14 G	10 cm
EBN-0114016	20 mm x 0.8 mm	14 G	16 cm
EBN-0114020	20 mm x 0.8 mm	14 G	20 cm
EBN-0116010	20 mm x 0.8 mm	16 G	10 cm
EBN-0116016	20 mm x 0.8 mm	16 G	16 cm
EBN-0116020	20 mm x 0.8 mm	16 G	20 cm
EBN-0118010	20 mm x 0.8 mm	18 G	10 cm
EBN-0118016	20 mm x 0.8 mm	18 G	16 cm
EBN-0118020	20 mm x 0.8 mm	18 G	20 cm

EXPANDER[®]

ESOPHAGEAL DILATOR SET

Product Features

- Set of Eight & Sixteen Dilators
- Total Length 85 cm
- Marking Every 1 cm
- Built-in Radiopaque Metal Ring
- Highly Smooth and Flexible





Cat No.	Set Includes	Dilator Sizes
EDS-0108085	Eight Dilators + One Spring Tip Wire	5,7,9,11,13,15, 17,19 mm x 85 cm
EDS-0116085	Sixteen Dilators + Spring Tip Wire	5,6,7,8,9,10,11,12,13,14,15,16,17,18, 19,20 mm x 85 cm

PALLIATOR® - ESOPHAGEAL STENT

Product Features

- Nitinol Stent Highly Radiopaque
- Covered From Inside
- Easy Retrieval Using Drawstring
- Dynamic Radial Force for Maximum Lumen Size
- Dumbbell Shape for Anti-Migration
- Flexible Construction to Conform to Anatomy
- X-Ray Visible Marks
- Scaled Marking on Delivery System




Cat No.	Stent OD	Flare OD	Stent Length	Delivery Size
EMS-0118008	18 mm	22 mm	8 cm	7 mm x 75 cm
EMS-0118010	18 mm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012	18 mm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015	18 mm	22 mm	15 cm	7 mm x 75 cm
EMS-0118010(AR)	18 mm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012(AR)	18 mm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015(AR)	18 mm	22 mm	15 cm	7 mm x 75 cm

*Other sizes are available on request

PALLIATOR[®] - BILIARY STENTS

Product Features

- Nitinol Stent
- Small Cell Size
- Dynamic Radial Force
- Braided Tubing of Delivery System
- Flexible Construction to Confirm Anatomy
- Zero Shortening at Proximal end
- Three Marking on Delivery System for Reliable Positioning





Cat No.	Stent Type	Stent OD	Stent Length	Delivery Size
EMS-0210006	Uncovered	10 mm	6 cm	8Fr x 180cm
EMS-0210007	Uncovered	10 mm	7 cm	8Fr x 180cm
EMS-0210008	Uncovered	10 mm	8 cm	8Fr x 180cm
EMS-0210009	Uncovered	10 mm	9 cm	8Fr x 180cm
EMS-0210010	Uncovered	10 mm	10 cm	8Fr x 180cm
EMS-0210011	Uncovered	10 mm	11 cm	8Fr x 180cm
EMS-0210012	Uncovered	10 mm	12 cm	8Fr x 180cm
EMS-0310006	Covered	10 mm	6 cm	10Fr x 180cm
EMS-0310007	Covered	10 mm	7 cm	10Fr x 180cm
EMS-0310008	Covered	10 mm	8 cm	10Fr x 180cm

PALLIATOR[®] Pseudocyst Stent

Product Features

- Nitinol Stent
- Fully Covered
- With Lasso
- Braided Tubing of Delivery System
- Three Marking on Delivery System for Reliable Positioning




Cat No.	Stent Type	Stent OD	Stent Length	Delivery Size
EMS-0414020	Covered	14 mm	2 cm	10Fr x 180cm
EMS-0414030	Covered	14 mm	3 cm	10Fr x 180cm
EMS-0414040	Covered	14 mm	4 cm	10Fr x 180cm
EMS-0416020	Covered	16 mm	2 cm	10Fr x 180cm
EMS-0416030	Covered	16 mm	3 cm	10Fr x 180cm
EMS-0416040	Covered	16 mm	4 cm	10Fr x 180cm

PALLIATOR[®] - DUODENAL / PYLORIC / COLONIC STENT

Product Features

- Nitinol Stent
- Small Cell Size
- Zero Shortening at Proximal End
- Braided Tubing of Delivery System
- Flexible Construction to Confirm Anatomy
- With Lasso
- Three Marking on Delivery System for Reliable Positioning



Cat No.	Stent Type	Stent OD	Stent Length	Delivery System
EMS-0522008	Uncovered	22 mm	8 cm	10Fr x 220cm
EMS-0522010	Uncovered	22 mm	10 cm	10Fr x 220cm
EMS-0522012	Uncovered	22 mm	12 cm	10Fr x 220cm
EMS-0524008	Uncovered	24 mm	8 cm	10Fr x 220cm
EMS-0524010	Uncovered	24 mm	10 cm	10Fr x 220cm
EMS-0524012	Uncovered	24 mm	12 cm	10Fr x 220cm
EMS -0622008	Covered	22 mm	8 cm	10Fr x 220 cm
EMS -0622010	Covered	22 mm	10 cm	10Fr x 220 cm
EMS-0622012	Covered	22 mm	12 cm	10Fr x 220 cm
EMS-0624008	Covered	24 mm	8 cm	10Fr x 220 cm
EMS-0624010	Covered	24 mm	10 cm	10Fr x 220 cm
EMS-0624012	Covered	24 mm	12 cm	10Fr x 220 cm



Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy

GI Endoscopy
GI Endoscopy


GI Endoscopy

Track-Wire[®]

GUIDE-WIRE

Product Features

- Nitinol Core
- Radiopaque Distal-Tungsten Filled Tip (5 cm)
- PTFE Jacket
- Hydrophilic Coated
- Twin Color




Cat No.	Diameter	Length
EGW- 0225450	0.025"	450 cm
EGW- 0232450	0.032"	450 cm
EGW- 0235450	0.035"	450 cm

*Other sizes are available on request

BILIARY DILATOR

Product Features

- Tapered and Rounded Ends
- Well Balanced Wall Thickness
- Fully Radiopaque




Cat No.	Tip OD	Dilator OD	Dilator Length
EBD- 0107200	4 Fr	7 Fr	200 cm
EBD- 0185200	5.5 Fr	8.5 Fr	200 cm
EBD- 0110200	6 Fr	10 Fr	200 cm

Smart-CutTM

SPHINCTEROTOMES

Product Features

- Monofilament Cutting Wire
- Triple Lumen Version
- Color Marking
- Separate Port for Contrast & Guidewire
- Pre-curved
- Ergonomical Shape Handle




Cat No.	OD	Cutting wire	Type
ESP- 0124020	2.4 mm	20 mm	Standard
ESP- 0124030	2.4 mm	30 mm	Standard
ESP- 0224007	2.4 mm	7 mm	Needle Knife
ESP- 0324020	2.4 mm	20 mm	Insulated Wire
ESP- 0324030	2.4 mm	30 mm	Insulated Wire

Drain-XTM

DRAINAGE CATHETER - INTRODUCERS - NBD CATHETER

Product Features

- Various Type and Sizes Available
- Tapered and Rounded Ends
- Well Balanced Wall Thickness
- Radiopaque Material
- Flaps or Pigtails
- Single Step 8.5 Fr & 10 Fr Introducer





Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy

GI Endoscopy
GI Endoscopy

GI Endoscopy

Catalog No.					Size
Biliary Amsterdam/ Curved Type	Biliary Double Pigtail	Pancreatic Curved	Biliary / Pancreatic Single Pigtail	Biliary Multiflap/Tannenbaum	
		EBS - 0305005	EBS - 0405005		5 Fr X 5 cm
		EBS - 0305007	EBS - 0405007		5 Fr X 7 cm
		EBS - 0305010	EBS - 0405010		5 F X 10 cm
EBS - 0107005	EBS - 0207005	EBS - 0307005	EBS - 0407005		7 Fr X 5 cm
EBS - 0107007	EBS - 0207007	EBS - 0307007	EBS - 0407007		7 Fr X 7 cm
EBS - 0107010	EBS - 0207010	EBS - 0307010	EBS - 0407010		7 Fr X 10 cm
EBS - 0107012	EBS - 0207012	EBS - 0307012	EBS - 0407012		7 Fr X 12 cm
EBS - 0107015	EBS - 0207015		EBS - 0407015		7 Fr X 15 cm
EBS - 0185005	EBS - 0285005		EBS - 0485005	EBS - 0585005	8.5 Fr x 5 cm
EBS - 0185007	EBS - 0285007		EBS - 0485007	EBS - 0585007	8.5 Fr x 7 cm
EBS - 0185010	EBS - 0285010		EBS - 0485010	EBS - 0585010	8.5 Fr x 10 cm
EBS - 0185012	EBS - 0285012		EBS - 0485012	EBS - 0585012	8.5 Fr x 12 cm
EBS - 0185015	EBS - 0285015		EBS - 0485015	EBS - 0585015	8.5 Fr x 15 cm
EBS - 0110005	EBS - 0210005		EBS - 0410005	EBS - 0510005	10 Fr X 5 cm
EBS - 0110007	EBS - 0210007		EBS - 0410007	EBS - 0510007	10 Fr X 7 cm
EBS - 0110010	EBS - 0210010		EBS - 0410010	EBS - 0510010	10 Fr X 10 cm
EBS - 0110012	EBS - 0210012		EBS - 0410012	EBS - 0510012	10 Fr X 12 cm
EBS - 0110015	EBS - 0210015		EBS - 0410015	EBS - 0510015	10 Fr X 15 cm
ENB - 0107250	NBD Catheter Pigtail with loop				7 Fr X 250 cm
EBS - 1005220	Pusher				5 Fr X 220 cm
EBS - 1007220	Pusher				7 Fr X 220 cm
EBS - 1085220	Introducer Set				8.5 Fr X 220 cm
EBS - 1010220	Introducer Set				10 Fr X 220 cm

*Other length e.g. 4 cm, 6cm, 8cm, 9cm, 11cm, 13cm can be made on order

**All sizes are available in sets also. Please add "S" after Catalogue Nos "e.g. : EBS - XXXXXXX-S"



Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy


GI Endoscopy
GI Endoscopy

GI Endoscopy

STONE EXTRACTION BASKET

Product Features

- Four Wire Construction
- Braided Wires of Hard Steel
- Separate Flushing / Contrast Port



Cat No.	Sheath OD / Length	Basket OD	No. of Wires
EBS - 0218020	1.8 mm / 230 cm	20 mm	Four Wire
EBS - 0224025	2.4 mm / 230 cm	25 mm	Four Wire
EBS - 0224030	2.4 mm / 230 cm	30 mm	Four Wire



STONE EXTRACTION BALLOON

Product Features

- Triple Lumen Variant
- Multi Diameter Balloon Inflation
- X-Ray Visible Marks
- Very Strong Balloon Material



Cat No.	Lumen	Sheath OD	Injection	Inflated Diameter
EEB- 0123200	Triple Lumen	2.8 mm	Above	11.5mm - 16mm
EEB- 0223200	Triple Lumen	2.8 mm	Below	11.5mm - 16mm



Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy

GI Endoscopy


GI Endoscopy

GI Endoscopy

LITHOTRIPTER - MANUAL

Product Features

- Light Weight
- Easy Rotation
- Compatable with all Types of Basket
- Good for Impacted Stone



Cat No.	Discription	Sheath OD	Length
EML-0160090	Handle with Spiral	6 mm	90 cm
EML-0325200	4 Wire Basket	2.3 mm	200 cm
EML-0330200	4 Wire Basket	2.3 mm	200 cm
EML-0560090	Metal Spiral	6 mm	90 cm

LITHOTRIPTER - ENDOSCOPIC

Product Features

- Easy to Use
- Light Weight
- All Component Aavailable Separately
- Through the Scope Working Channel



Cat No.	Discription	OD	Length
EML-0237210	Handle + Spiral + Basket	3.7 mm	210 cm
EML-0425210	4 Wire Basket	25 mm	210 cm
EML-0430210	4 Wire Basket	30 mm	210 cm
EML-0537210	Metal Spiral	3.7 mm	210 cm
EML-0625210	Metal Spiral with Basket	3.7 mm / 25 mm	210 cm
EML-0630210	Metal Spiral with Basket	3.7 mm / 30 mm	210 cm



Endo-Med[®]
A Resurgence to life

GI Endoscopy

GI Endoscopy


GI Endoscopy
GI Endoscopy

GI Endoscopy

DISPOSABLE CYTOLOGY BRUSH

Product Features

- Nitinol Core
- Radiopaque Distal-Tungsten Filled Tip (5 cm)
- PTFE Jacket
- Hydrophilic Coated
- Twin Color



Cat No.	Diameter	Length
ECY - 0123180	2.3 mm	180 cm
ECY - 0123230	2.3 mm	230 cm

*Other sizes are available on request

SCOPE CLEANING BRUSH

Product Features

- Tapered and Rounded Ends
- Well Balanced Wall Thickness
- Fully Radiopaque




Cat No.	For Chanal Size	Brush	Length
ECB - R123180	2.8 mm - 4.2 mm	One Side	180 cm
ECB - R223230	2.8 mm - 4.2 mm	Both Side	230 cm

P.E.G. Set

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY FEEDING SET

Product Features

- Pull Technique
- Thin Wall Tubing
- Large Inner Lumen




Cat No.	OD	Kit Contents
EPG- 0114030	14 fr	Mini Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire
EPG- 0120030	20 fr	Mini Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire
EPG- 0124030	24 fr	Mini Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire
EPG- 0214030	14 fr	Full Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire, Snare, Scissors, Hemostat, Gauze, Drape, Needle, Syringe, Povidon Iodine,
EPG- 0220030	20 fr	Full Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire, Snare, Scissors, Hemostat, Gauze, Drape, Needle, Syringe, Povidon Iodine,
EPG- 0224030	24 fr	Full Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire, Snare, Scissors, Hemostat, Gauze, Drape, Needle, Syringe, Povidon Iodine,


Gastrostomy Feeding Tube

Product Features


- Medical Grade 100% Silicon Tube
- Thin Walls Tube
- Large Inner Lumen
- Non-traumatic Distal Tip



Cat No.	Tube OD	Tube Length	Balloon Volume
EPG-2112024	12 Fr	24 cm	5 ml
EPG-2114024	14 Fr	24 cm	5 ml
EPG-2116024	16 Fr	24 cm	20 ml
EPG-2118024	18 Fr	24 cm	20 ml
EPG-2120024	20 Fr	24 cm	20 ml
EPG-2122024	22 Fr	24 cm	20 ml
EPG-2124024	24 Fr	24 cm	20 ml



Pneumatic Inflation Device



Channel Cap with Flushing Tube



Spring Tip Wire



Balloon Inflation Device



Endoscope Bite Guard



 Endo-Med Technologies Pvt. Ltd.

Corporate Off. : Off. No. 607, 6th Floor, DDA Building,
Plot No. 4, Laxmi Nagar Distt. Centre, Delhi-110092. INDIA

Factory : A-138, Sector 63, Noida - 201301 U.P.

Phone: +91-11-22048222, 42448449

Customer Care No. :  +91-7840070477

E-mail: sales@endomed.in, **Web :** www.endomed.in

 Endo-Med Technologies Pvt. Ltd.  endomed_technologies

Distributed By :