


An ISO 13485:2016 and
ISO 9001:2015 Certified Company




Laparoscopy

ELECTRO SURGICAL INSTRUMENT

Product Features

- Ergonomically designed soft handle for comfort while using
- Easy one figure rotatable
- Very light weight
- Easy operation of ratchet




Cat No.	Specification	Size
ELH-0105033	Scissor Straight Single Use	5 mm X 330 mm
ELH-0205033	Scissor Metzenbaum Curved Single Use	5 mm X 330 mm
ELH-0305033	Scissor Fully Insulated Single Use	5 mm X 330 mm
ELH-0405033	Dissector Maryland Single Use	5 mm X 330 mm
ELH-0505033	Dissector Mixer Single Use	5 mm X 330 mm
ELH-0605033	Dissector Straight Single Use	5 mm X 330 mm
ELH-0705033	Grasper Fenestrated Double Action Single Use	5 mm X 330 mm
ELH-0805033	Grasper Fenestrated Single Action Single Use	5 mm X 330 mm
ELH-0905033	Grasper Clinch Single Use	5 mm X 330 mm
ELH-1005033	Grasper Babcock Single Use	5 mm X 330 mm
ELH-1105033	Grasper Allis Single Use	5 mm X 330 mm
ELH-1205033	Grasper Rat-Teeth Single Use	5 mm X 330 mm
ELH-1305033	Grasper Blunt Jaw Single Use	5 mm X 330 mm

LAPROSCOPIC TROCAR

Product Features

- Single use only, avoid cross contamination
- Unique design, small trauma, quick recovery
- Thread design, ovoids movement of trocar sleeve during operation.



Cat No.	Specification	Size
ELT-S103100	Trocar Optical Single Use	3 mm x 100 mm
ELT-S103150	Trocar Optical Single Use	3 mm x 150 mm
ELT-S105100	Trocar Optical Single Use	5 mm x 100 mm
ELT-S105150	Trocar Optical Single Use	5 mm x 150 mm
ELT-S110100	Trocar Optical Single Use	10 mm x 100 mm
ELT-S110150	Trocar Optical Single Use	10 mm x 150 mm
ELT-S112100	Trocar Optical Single Use	12 mm x 100 mm
ELT-S112150	Trocar Optical Single Use	12 mm x 150 mm
ELT-S205100	Trocar Bladed Single Use	5 mm x 100 mm
ELT-S205150	Trocar Bladed Single Use	5 mm x 150 mm
ELT-S210100	Trocar Bladed Single Use	10 mm x 100 mm
ELT-S210150	Trocar Bladed Single Use	10 mm x 150 mm
ELT-S212100	Trocar Bladed Single Use	12 mm x 100 mm
ELT-S212150	Trocar Bladed Single Use	12 mm x 150 mm

RETRIEVAL BAG

Product Features

- The bag body of high strength with good biocompatibility.
- TPU material, Transparent soft bag body
- Variety of bag sizes to choose from
- Automatic tightening of mouth, to prevent overflow




Cat No.	Specification	Size
ERB-0110250	Retrieval Bag Single Use	250 ml
ERB-0110500	Retrieval Bag Single Use	500 ml
ERB-0110750	Retrieval Bag Single Use	750 ml
ERB-0110999	Retrieval Bag Single Use	1000 ml
ERB-0210250	Retrieval Bag Single Use (Manual Type)	250 ml
ERB-0210500	Retrieval Bag Single Use (Manual Type)	500 ml
ERB-0210750	Retrieval Bag Single Use (Manual Type)	750 ml
ERB-0210999	Retrieval Bag Single Use (Manual Type)	1000 ml

SUCTION & IRRIGATION SET

Product Features

- Ergonomically designed handpiece can be easily operated by either hand
- Stainless steel probe with position function, easy to fix and draw out
- Flip designed rear port for easy access inside chamber
- Transparent headpiece, unobstructed visualization
- Clear channel, consistent diameter with self-flush function



Cat No.	Specification	Size
ESI-0105360	Suction & Irrigation Set Trumpet Type	5 mm x 360 cm
ESI-0105430	Suction & Irrigation Set Trumpet Type	5 mm x 430 cm

DISPOSABLE VERESS NEEDLE

Product Features

- Stainless steel needle with strong polycarbonate handle
- Double warning system confirms proper insertion
- Ergonomic handle for easy controlled insertion



Cat No.	Specification	Size
EVN-0121120	Veress Needle Single Use OD	2.1mm X 120mm
EVN-0121150	Veress Needle Single Use OD	2.1mm X 150mm

LAPAROSCOPY ELECTRODES

Product Features

- Monopolar Electrode
- Insulated Tip



Cat No.	Items	Size
ELE-0105330	J Hook Electrode	5mm X 330 mm
ELE-0205530	L Hook Electrode	5mm X 330 mm
ELE-0305330	Spatula Electrode	5mm X 330 mm
ELE-0405330	Needle Electrode	5mm X 330 mm



Laparoscopy Laparoscopy Laparoscopy

Company Profile

Endo-Med Technologies Pvt. Ltd. is having its corporate office at New Delhi, capital city of India and manufacturing facility at Noida (National Capital Region) engaged in development, manufacturing and globally distribution of minimal invasive medical devices in the field of Gastro Endoscopy, Laparoscopy, and Bronchoscopy. Backed by more than 25 years of experience, Endo-Med is committed to serve the customers with quality products at affordable price and service to the utmost satisfaction of the customers.

Quality Policy


The members of "team Endo-Med" are experienced with the largest companies across the world and exposed to trans-continent markets, thus possessing the highest level of expertise in their respective fields. The company is an ISO 13485:2016 certified organization which is deeply involved in the manufacturing of the best quality medical products under the best quality management systems. With guidance from the renowned Doctors in India and abroad associated with Endo-Med, we develop the product keeping in mind, the every minute expectation of the user from the product thus making our product among the most user friendly for a medical practitioner in their procedural room. All Endo-Med Products are "CDSCO" Certified and are manufactured as per the guidelines of GMP*.

Philosophy

At Endo-Med, Quality is the way of life and operations. Stage wise rigorous inspection is carried out right from procurement of raw material to delivery of the finished product. All of our employees are committed to constantly improve the quality management system.

The company believes in long term relationship and mutual growth as objective. At Endo-Med, quality is a culture and service a tradition.

Endo-Med Manufacturing Facility at Noida



Endo-Med Technologies Pvt. Ltd.

Corporate Off. : Off. No. 607, 6th Floor, DDA Building, Plot No. 4, Laxmi Nagar Distt. Centre, Delhi-110092. INDIA

Factory : A-138, Sector 63, Noida - 201301 U.P.

Phone: +91-11-22048222, 42448449

Customer Care No. : +91-7840070477

E-mail: sales@endomed.in, **Web :** www.endomed.in

Facebook: Endo-Med Technologies Pvt. Ltd. **Instagram:** endomed_technologies

Distributed By :